

## DH OFFLOADING POST-OP STYLE SHOE

Designed for use in the treatment of plantar wounds, this shoe can be an excellent alternative to post-op style shoes with static pads (wool, felt, foam, etc.). It accommodates the ambulatory and compliant patient for short-term treatment (1 to 4 months).

### INDICATIONS FOR USE

- Plantar wounds
- Post-operative use
- Plantar foot ulcers

### FEATURES AND BENEFITS

- Designed to accommodate either left or right foot
- Designed around Össur's unique Active Offloading Insole
- Designed to accommodate bulky dressings
- Super soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe pieces aides in wound accessibility for daily wound site inspection



### REQUIRED MEASUREMENTS

- Shoe size

## ORDERING INFORMATION

	Shoe Size UK		Shoe Size EU		
Size	Male	Female	Male	Female	Part#
Small	3.5 - 5.5	4 - 6	35.5 - 38	37 - 39	10340
Medium	6 - 8	6.5 - 8.5	38.5 - 42	39.5 - 43.5	10341
Large	8.5 - 10.5	9 - 10	42.5 - 45	44 - 45	10342
X-Large	11 - 13.5	10.5+	45.5 - 48.5	45.5-48.5	10343

## ACCESSORIES

### REPLACEMENT PARTS

Size	Insoles
Small	12365
Medium	12366
Large	12367
X-Large	12368