

Footwear

For Healthy Walking

FREE XL Series

Specialized Footwear for People with Diabetes and Peripheral Neuropathy

Key Features & Benefits:

- Seamless joint area to minimize irritation
- Optimized rocker angle (~20°) for improved foot movement efficiency
- non-slip rubber track
- Includes two adjustment inlays to regulate interior volume
- Carbon spring reinforcement over partial length of the foot plate.
- 8mm wider across the metatarsal heads, plus 2 x 5mm EVA removable footbeds.

MS SC010

0200 Black

Fully opening sandal.
Carbon Spring System
with gait control.

Free XL

35-46

2x5 mm EVA

KG

≈386 g/
Unit (size 38)

Artificial
Leather

Textile

EVA/Rubber
w/ Carbon Insole

Well-padded ankle area
for enhanced comfort

Turn-return Velcro straps
with adjustable system
for easy on-off



Read the [scientific paper](#)
supporting the Free XL design

Rehabilitation Footwear

For Healthy Walking

FREE XL Series

MS SC009

0200 Black

Wide opening shoe.
Carbon Spring System
with gait control.



Soft and padded
textile tongue

Well-padded ankle area
for enhanced comfort



Turn-return Velcro straps
with adjustable system
for easy on-off

Loops for easy
donning and doffing

For more information, orders, or inquiries:

Orthotic Prosthetic Supplies Ltd

Unit 73, Basepoint Business Centre
Metcalf Way, Crawley
RH11 7XX

Tel: 07966656977

Email: info@opsupplies.co.uk